

## MHS Hospital Bed Weighing Scales

MHS are high precision digital hospital bed scales which allow bedridden patients to be simply and accurately weighed with no or little disturbance.

These bed scales are used primarily in hospitals for weighing acute stroke and other patients.



Portable Trolley Bed Scale



Weighing Pad

The main features of the MHS scales are:

- Suits any bed size
- Portable trolley design
- Powered by rechargeable battery and AC/DC power adapter
- Simple to setup and use, little effort required to roll bed casters on the 4 weighing pads
- Weighing units: kilograms (kg) and pounds (lb)
- High capacity 500kg x 0.1kg / 1100lb x 0.2lb
- Weight of the each bed type could be stored in memory and recalled when needed
- Tare (weigh of the bed) could be entered via keypad before or when the bed is on the scale
- Light weight 4 weighing pads
- Easy to read large 0.8"/ 20 mm 7-segment red LED display
- LED annunciators provide indication of the weighing and tare status
- Optional RS232 serial data communication port
- Made in Canada
- Warranty: 1 year

Mars Scale Manufacturing Inc.

697 Petrolia Road, North York, Ontario, Canada M3J 2N6

Phone: 416-667-9557 Fax: 416-667-9943

Website: [www.marsscale.com](http://www.marsscale.com) E-mail: [info@marsscale.com](mailto:info@marsscale.com)